



CELLECTA
USER MANUAL

CRISPRuTest™

Functional Cas9

Activity Assay Kit

v1 — Last update: 2019/01/10

Cellecta, Inc.

Table of Contents

- 1. CRISPRuTest™ Functional Cas9 Activity Assay Kit for Any Cell (Mammalian)..... 1
- 2. Kit Components 2
- 3. Additional Required Materials 3
- 4. Protocol Overview 4
 - 4.1. Assay Procedure 4
 - 4.2. Cas9 Activity Calculation 7
- 5. Technical Support 8
- 6. Safety Guidelines 9
- 7. Contact Us 10

1. CRISPRuTest™ Functional Cas9 Activity Assay Kit for Any Cell (Mammalian)

The CRISPRuTest™ Functional Cas9 Activity Assay Kit for Any Cell (Mammalian) (Cat. #CRUTEST) can be used to measure the Cas9 nuclease activity in any mammalian cell system. The CRISPR system uses a target single-guide RNA molecule (sgRNA) to target and cut a specific sequence on the genomic DNA with the Cas9 nuclease. The CRISPRuTest Kit provides a FACS-based assay to assess the Cas9 nuclease activity in cells expressing *Streptococcus pyogenes* Cas9. The kit contains pseudoviral packaged lentiviral constructs with a green fluorescent protein (GFP) gene whose fluorescence is disrupted in cells with active Cas9. The kit provides enough reagents to measure this activity in 5-15 cell lines.

References and Product Citations for all Cellecta products can be found on the Cellecta website:

<https://www.cellecta.com/resources/publications/>.

Please read the entire user manual before proceeding with your experiment. Also, please note that, when working with pseudoviral particles, you should follow the recommended guidelines for working with Biosafety Level 2 (BSL-2) materials.



Click the ? Download as PDF link located at the bottom of the left menu to download the PDF version of this user manual.

PLEASE NOTE: The purchase of all Cellecta Products is covered by Cellecta's standard Terms and Conditions of Sale as described on [Cellecta's website](#), and selected Products containing particular technology or having certain features are also subject to restrictions on use as outlined in the [Label License section](#) of our website. Please review these Terms and Label License Restrictions before opening and using your Product and, if you are not able to abide by the restriction, contact Cellecta to return the item to Cellecta for a full refund.

Last modified: 2019/01/10

2. Kit Components

Component	Volume
CRISPRuTest™ CT-Active (CT-A) Viral Mix	500 µl
CRISPRuTest™ CT-Background (CT-B) Viral Mix	500 µl
Transduction Reagent (1000X)	50 µl

The CRISPRuTest™ Kit should be stored at -80°C until ready for use.

Last modified: 2017/05/30

3. Additional Required Materials

Use of the CRISPRuTest assay requires access to a flow cytometer with the following excitation and emission requirements:

Excitation: 488nm and 561nm

Emission: 530/20nm and 590/20nm

Cellecta also offers the following products that may be useful when running this assay:

Cas9-Only Plasmids

Cat.#	Description	Quantity	Price
SVC9B-PS	CRISPR Cas9 pRCCB-CMV-Cas9-2A-Blast (plasmid)	25 µg	\$300
SVC9-PS	CRISPR Cas9 pRCCH-CMV-Cas9-2A-Hygro (plasmid)	25 µg	\$300
SVC9G-PS	CRISPR Cas9 pRCCG-CMV-Cas9-2A-TagGFP2 (plasmid)	25 µg	\$300
SVC9R-PS	CRISPR Cas9 pRCCR-CMV-Cas9-2A-TagRFP (plasmid)	25 µg	\$300
SVC9P-PS	CRISPR Cas9 pRCCP-CMV-Cas9-2A-Puro (plasmid)	25 µg	\$300
SVC9B-VS	CRISPR Cas9 pRCCB-CMV-Cas9-2A-Blast (virus)	1 × 10 ⁶ TU	\$500
SVC9-VS	CRISPR Cas9 pRCCH-CMV-Cas9-2A-Hygro (virus)	1 × 10 ⁶ TU	\$500
SVC9G-VS	CRISPR Cas9 pRCCG-CMV-Cas9-2A-TagGFP2 (virus)	1 × 10 ⁶ TU	\$500
SVC9R-VS	CRISPR Cas9 pRCCR-CMV-Cas9-2A-TagRFP (virus)	1 × 10 ⁶ TU	\$500
SVC9P-VS	CRISPR Cas9 pRCCP-CMV-Cas9-2A-Puro (virus)	1 × 10 ⁶ TU	\$500

Positive Control Cell Line

Cas9+ MDA-MB-231 Cells may be run in parallel with your target cells as a positive control for the assay.

Cat.#	Description	Quantity	Price
ZMDAMB231-C9	Cas9+ MDA-MB-231 Cells with Blast selection	1 × 10 ⁶ cells	\$1,450

Transduction Control Virus

Non-Targeting CRISPR Control pRSG16-U6-(sg)-UbiC-TagRFP-2A-Puro (virus) may be used to assess transduction efficiency of the target cells.


Cat.#	Description	Quantity	Price
SGCTL-NT-pRSG16-V	Non-Targeting CRISPR control pRSG16-U6-(sg)-UbiC-TagRFP-Puro Vector (virus)	1 × 10 ⁷ TU	\$500

Other than the specific reagents and instruments listed above, the protocols assume the user has access to standard materials (e.g., polypropylene tubes, pipette tips), equipment (table top centrifuges, pipettes, scales), and common reagents (e.g., TE buffer, ethanol) and buffers used in a typical life science laboratory.

Last modified: 2018/07/04

4. Protocol Overview

Both the CT-Active (CT-A) and CT-Background (CT-B) reagents contain a lentivirus expressing RFP and GFP from one transcript, separated by a T2A peptide linker sequence (see Figure below). Transduction of either of these reagents into cells will initially produce cells expressing both GFP and RFP fluorescence.



CRISPRuTest Vector Diagrams CT-A and CT-B

The difference between the two viral reagents is that the CT-A expresses a GFP-targeting sgRNA, whereas the CT-B viral reagent expresses a non-targeting sgRNA. Therefore, within a few days after transduction, CT-A transduced cells will lose GFP and only have RFP fluorescence, whereas CT-B transduced cells will continue to exhibit both RFP and GFP fluorescence. The decline in GFP fluorescence in the CT-A cell population, then, can be used to assess Cas9 activity and forms the basis of the CRISPRuTest assay.

Last modified: 2017/05/23

4.1. Assay Procedure

Two populations of putative Cas9-positive cells are transduced with the CRISPRuTest™ CT-A and CT-B viral mix reagents, respectively. After 4 days, each population is analyzed by flow by cytometry. The disappearance of GFP in the CT-A population as compared with the CT-B population provides a quantitative measure of the Cas9 activity in the target cells.

The assay was optimized using MDA-MB-231 cells. These Cas9-expressing MDA-MB-231 with blasticidin selection marker cells are available from Cellecta for use in the assay as a positive control (see [Additional Required Materials](#)). Some optimization may be needed based on the growth characteristics of your target cells. If a negative control is also desired, include a cell line not expressing Cas9 (ideally, the parental cells for the Cas9-positive target cells) in a parallel run of the assay.

Day 0

1. Quickly thaw the CRISPRuTest™ CT-A and CT-B lentiviral particles in a water bath at 37°C. Transfer the thawed particles to a laminar flow hood, gently mix by rotation, inversion, or gentle vortexing, and keep on ice. Unused reagent can be aliquoted, refrozen at -80°C, and used again for subsequent experiments.

2. Suspend Cas9+ cells in growth medium supplemented with 1X Transduction Reagent, at a density of ca. 100,000 cells/ml.

Note: This cell density was calculated for MDA-MB-231 cells. Depending on cell size and growth, you may need to use a different concentration and correspondingly-sized plate. As a rule of thumb, cells should be transduced at a density such that they would become confluent in ~48 hours. For the assay, you should plate at least 100,000 cells.

3. Aliquot 1 ml of cell suspension (100,000 cells) into each of 2 wells of a 12-well plate.
4. Add 20 µl of CT-A virus into one well and 20 µl of CT-B virus into the other well, and then mix and return cells to incubator.

Note: For most cell lines, 20 µl of CT viral reagents will suffice to obtain at least 50% RFP+ cells (the recommended minimum percentage of transduced cells for optimal assay sensitivity). For hard-to-transduce cell lines, more virus might be needed. If in doubt, it is recommended to set up two sets of transductions with 20 µl of CT-A and CT-B as described above, and 50 µl of each in the second transductions. For the final calculation, use the samples that have 50%-80% RFP+ cells on Day 4. A control lentiviral vector may be used before running the assay to test the transduction efficiency of your cells (see [Additional Required Materials](#)).

Days 1-5

Exchange medium with fresh growth medium and grow cells under standard conditions for 3 days. Passage cells as needed. Cells should not become confluent.

Day 6


1. Analyze cells by flow cytometry, using settings below:

Channel 1: excitation 488nm, emission 530/20nm (GFP)


Channel 2: excitation 561nm, emission 590/20nm (RFP)

Note: Cells will typically show a significant loss of GFP activity as soon as Day 4, but cells with lower Cas9 levels may only show marginal effects in this time period. For more robust data, we recommend incubating cells for 6 days.


2. For both CT-B and CT-A samples, gate in RFP+ cells as shown below:



#_Plot gated cells of CT-B sample GFP vs. RFP as shown below, adjusting Channel 1 and 2 intensity so that GFP and RFP signal falls into a diagonal. Then, gate all cells below diagonal (“GFP_low” cells):



3. Do the same for CT-A sample, applying the same “GFP_low” gate:



Last modified: 2017/11/10

4.2. Cas9 Activity Calculation

The difference between the percentage of “GFP_low” cells in the CT-A population and percentage of “GFP_low” cells in the CT-B population indicates the percentage of cells which underwent Cas9-mediated GFP knockout in the CT-A sample. The percentage of “GFP_low” cells in the CT-B population should be close to 0% (i.e., no sgRNA targeting GFP), where as in the CT-A population it should be near 100% (sgRNA targeting GFP).

The flow cytometry data from the the previous section using the positive control cell line MDA-MB-231-Cas9-Blast (see [Additional Required Materials](#) to obtain this cell line) can be used to calculate the percentage of GFP knockout, which is a direct measurement of the Cas9 activity, as follows:

Target Cell Population	“GFP low”
CT-A Cells	0.938
CT-B Cells	0.0162

Percentage GFP Knockout in Target Cas9 cells = 93.8% – 1.62% = 92.2%

Note: For a Cas9-negative cell line, the percentage of “GFP_low” cells in both CT-A and CT-B should be similarly close to 0%.

Last modified: 2017/05/23

5. Technical Support

Email Addresses

Technical Support: tech@cellecta.com

General information: info@cellecta.com

Phone Numbers

Phone: +1 650 938-3910

Toll-free: (877) 938-3910

For the latest technical news and updates, visit Cellecta's blog at: <https://www.cellecta.com/blog-news/>

Last modified: 2017/05/08

6. Safety Guidelines

The HIV-based lentivector system is designed to maximize its biosafety features, which include:

- A deletion in the enhancer of the U3 region of 3'ΔLTR ensures self-inactivation of the lentiviral construct after transduction and integration into genomic DNA of the target cells.
- The RSV promoter upstream of 5'LTR in the lentivector allows efficient Tat-independent production of lentiviral RNA, reducing the number of genes from HIV-1 that are used in this system.
- Number of lentiviral genes necessary for packaging, replication and transduction is reduced to three (gag, pol, rev). The corresponding proteins are expressed from different plasmids lacking packaging signals and share no significant homology to any of the expression lentivectors, pVSV-G expression vector, or any other vector to prevent generation of recombinant replication-competent virus.
- None of the HIV-1 genes (gag, pol, rev) will be present in the packaged lentiviral genome, as they are expressed from packaging plasmids lacking packaging signal—therefore, the lentiviral particles generated are replication-incompetent.
- Lentiviral particles will carry only a copy of your expression construct.

Despite the above safety features, use of HIV-based vectors falls within NIH Biosafety Level 2 criteria. For a description of laboratory biosafety level criteria, consult the Centers for Disease Control Office of Health and Safety Web site at:

<https://www.cdc.gov/biosafety/publications/bmb15/index.htm>

It is also important to check with the health and safety guidelines at your institution regarding the use of lentiviruses and follow standard microbiological practices, which include:

- Wear gloves and lab coat at all times when conducting the procedure.
- Always work with lentiviral particles in a Class II laminar flow hood.
- All procedures are performed carefully to minimize the creation of splashes or aerosols.
- Work surfaces are decontaminated at least once a day and after any spill of viable material.
- All cultures, stocks, and other regulated wastes are decontaminated before disposal by an approved decontamination method such as autoclaving. Materials to be decontaminated outside of the immediate laboratory area are to be placed in a durable, leakproof, properly marked (biohazard, infectious waste) container and sealed for transportation from the laboratory.

Last modified: 2018/06/03

7. Contact Us

Email Addresses

General information: info@cellecta.com

Sales: sales@cellecta.com

Orders: orders@cellecta.com

Technical Support: tech@cellecta.com

Mailing Address

Cellecta, Inc.
320 Logue Ave.
Mountain View, CA 94043
USA

Phone Numbers

Phone: +1 650 938-3910

Toll-free: (877) 938-3910

Fax: +1 650 938-3911

For the latest technical news and updates, visit Cellecta's blog at:

<https://www.cellecta.com/blog-news/>



Copyright © 2019 Cellecta, Inc. All Rights Reserved.

Last modified: 2019/01/10