AdvanBlock TM-PF Blocking Solution

Non-protein blocking solution for fluorescent and chemiluminescent Western blots

AdvanBlock-PF is a protein-free blocking solution, optimized for use with the WesternBright MCF fluorescent Western blotting kit and also an excellent choice for chemiluminescent Western blot experiments. This fast-acting blocking solution stabilizes the fluorescence of the WesternBright MCF secondary antibody conjugates. AdvanBlock-PF can reduce background when used with primary antibodies that have a high degree of cross-reactivity with protein blockers such as BSA, casein, or milk protein. With low-quality primary antibodies that may require protein-based blocking agents, BSA or non-fat dry milk can be dissolved directly in the AdvanBlock-PF solution used to dilute the antibody. Provided as a 5X concentrate.



Advantages

- PROTEIN-FREE for lower background with some antibodies
- VERSATILE compatible with fluorescent and chemiluminescent detection
- OPTIMIZED stabilizes fluorescence of WesternBright MCF fluorescent secondary antibodies

Ordering Information

Catalog Number	Product		Size
R-03023-D20	AdvanBlock™-PF Blocking Solution		200 ml
Related Produc	cts		
K-12021-010	WesternBright [™] MCF fluorescent Western blotting kit (Contains: secondary antibodies, washing solution, blocking solution, pre-cut PVDF membranes and background quenching sheets)		10 blots
K-12042-D10	WesternBright™ Quantum Western Blotting HRP Substrate	For 1000 cm ² membrane	100 ml
K-12043-D10	WesternBright™ Sirius Western Blotting HRP Substrate	For 1000 cm ² membrane	100 ml
K-12045-D20	WesternBright™ ECL Western Blotting HRP Substrate	For 2000 cm ² membrane	200 ml
L-07014-100	LucentBlue™ X-ray Film	8 x 10" sheets	100 sheets
L-07032-005	Incubation Tray, Traditional Black Body with Clear Lid	6 x 9 cm	5/pack
R-05051-250	APC Goat-Anti-Rabbit IgG Antibody		250 µl
R-05052-250	RPE-Goat-Anti-Mouse IgG Antibody		250 µl
R-05071-500	Goat Anti-Mouse IgG (H+L), HRP Conjugate		500 µl
R-05072-500	Goat Anti-Rabbit IgG (H+L), HRP Conjugate		500 µl

Advansta Corporation

2140 Bering Drive | San Jose, CA 95131

Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com Product information: www.advansta.com/AdvanBlock-PF Copyright © 2011-2014 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. AdvanBlockTM, LucentBlueTM and WesternBrightTM are trademarks of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

